

Product datasheet for **GA114928**

Human EID2B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EID2B CRISPRa kit - CRISPR gene activation of human EP300 interacting inhibitor of differentiation 2B
Symbol:	EID2B
Components:	GA114928G1, EID2B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114928G2, EID2B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114928G3, EID2B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96D98</u>
Synonyms:	EID-2B; EID-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)