

Product datasheet for **GA114933**

Human JSRP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	JSRP1 CRISPRa kit - CRISPR gene activation of human junctional sarcoplasmic reticulum protein 1
Symbol:	JSRP1
UniProt ID:	Q96MG2
Synonyms:	JP-45; JP45
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)