

Product datasheet for **GA114934**

Human MOBKL2A (MOB3A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MOB3A CRISPRa kit - CRISPR gene activation of human MOB kinase activator 3A
Symbol:	MOB3A
Components:	GA114934G1, MOBKL2A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114934G2, MOBKL2A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114934G3, MOBKL2A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96BX8
Synonyms:	MOB-LAK; MOB1C; MOBKL2A; moblak
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)