

Product datasheet for **GA114966**

Human DNAI3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR63 CRISPRa kit – CRISPR gene activation of human WD repeat domain 63
Symbol:	DNAI3
Components:	GA114966G1, DNAI3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114966G2, DNAI3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114966G3, DNAI3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IWG1</u>
Synonyms:	DIC3; NYD-SP29
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

