

## Product datasheet for **GA114986**

### Human TYW3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TYW3 CRISPRa kit - CRISPR gene activation of human tRNA-yW synthesizing protein 3 homolog
Symbol:	TYW3
UniProt ID:	<a href="#">Q6IPR3</a>
Synonyms:	C1orf171
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)