

Product datasheet for **GA114987**

Human C1orf173 (ERICH3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ERICH3 CRISPRa kit - CRISPR gene activation of human glutamate rich 3
Symbol:	ERICH3
Components:	GA114987G1, C1orf173 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114987G2, C1orf173 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114987G3, C1orf173 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5RHP9</u>
Synonyms:	C1orf173
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)