

Product datasheet for **GA114988**

Human **LRR1Q3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRR1Q3 CRISPRa kit - CRISPR gene activation of human leucine rich repeats and IQ motif containing 3
Symbol:	LRR1Q3
UniProt ID:	<u>A6PVS8</u>
Synonyms:	LRR1Q3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)