

Product datasheet for **GA114989**

Human TPRG1L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TPRG1L CRISPRa kit - CRISPR gene activation of human tumor protein p63 regulated 1 like
Symbol:	TPRG1L
Components:	GA114989G1, TPRG1L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114989G2, TPRG1L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114989G3, TPRG1L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T0D9</u>
Synonyms:	FAM79A; h-mover; mover
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)