

Product datasheet for **GA115004**

Human DNAH14 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DNAH14 CRISPRa kit - CRISPR gene activation of human dynein axonemal heavy chain 14
Symbol:	DNAH14
Components:	GA115004G1, DNAH14 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115004G2, DNAH14 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115004G3, DNAH14 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q0VDD8</u>
Synonyms:	C1orf67; Dnahc14; HL-18; HL18
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)