

Product datasheet for **GA115007**

Human ZNF648 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF648 CRISPRa kit - CRISPR gene activation of human zinc finger protein 648
Symbol:	ZNF648
Components:	GA115007G1, ZNF648 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115007G2, ZNF648 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115007G3, ZNF648 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5T619
Synonyms:	OTTHUMP00000035419; zinc finger protein 648
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)