

Product datasheet for **GA115008**

Human TEDDM1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEDDM1 CRISPRa kit - CRISPR gene activation of human transmembrane epididymal protein 1
Symbol:	TEDDM1
Components:	GA115008G1, TEDDM1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115008G2, TEDDM1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115008G3, TEDDM1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T9Z0</u>
Synonyms:	EDDM9; Epdd1; HE9; HEL-S-45e; TMEM45C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)