

Product datasheet for **GA115019**

Human KIST (UHMK1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	UHMK1 CRISPRa kit - CRISPR gene activation of human U2AF homology motif kinase 1
Symbol:	UHMK1
UniProt ID:	Q8TAS1
Synonyms:	KIS; KIST; P-CIP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)