

Product datasheet for **GA115028**

Human TSACC activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSACC CRISPRa kit - CRISPR gene activation of human TSSK6 activating cochaperone
Symbol:	TSACC
Components:	GA115028G1, TSACC gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115028G2, TSACC gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115028G3, TSACC gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96A04</u>
Synonyms:	C1orf182; SIP; SSTK-IP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)