

## Product datasheet for **GA115029**

### Human IQGAP3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	IQGAP3 CRISPRa kit - CRISPR gene activation of human IQ motif containing GTPase activating protein 3
Symbol:	IQGAP3
UniProt ID:	<u><a href="#">Q86VI3</a></u>
Synonyms:	MGC1947; MGC10170; MGC10831
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)