


## Product datasheet for **GA115046**

### Human VSTM2L activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | VSTM2L CRISPRa kit - CRISPR gene activation of human V-set and transmembrane domain containing 2 like                                                                                                                                                                                                         |
| Symbol:       | VSTM2L                                                                                                                                                                                                                                                                                                        |
| Components:   | GA115046G1, VSTM2L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115046G2, VSTM2L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115046G3, VSTM2L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q96N03</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | C20orf102; dj1118M15.2                                                                                                                                                                                                                                                                                        |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)