

## Product datasheet for **GA115057**

### Human GPR73B (PROKR2) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PROKR2 CRISPRa kit - CRISPR gene activation of human prokineticin receptor 2
Symbol:	PROKR2
UniProt ID:	<a href="#">Q8NFJ6</a>
Synonyms:	dj680N4.3; GPR73b; GPR73L1; GPRg2; HH3; KAL3; PKR2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)