

Product datasheet for **GA115058**

Human SLX4IP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLX4IP CRISPRa kit - CRISPR gene activation of human SLX4 interacting protein
Symbol:	SLX4IP
Components:	GA115058G1, SLX4IP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115058G2, SLX4IP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115058G3, SLX4IP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5VYV7</u>
Synonyms:	bA204H22.1; bA254M13.1; C20orf94; dj1099D15.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)