


Product datasheet for **GA115065**

Human **SPLUNC3 (BPIFA3)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	BPIFA3 CRISPRa kit - CRISPR gene activation of human BPI fold containing family A member 3
Symbol:	BPIFA3
Components:	GA115065G1, SPLUNC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115065G2, SPLUNC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115065G3, SPLUNC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BQP9</u>
Synonyms:	MGC44525; SPLUNC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)