

Product datasheet for **GA115079**

Human TMEM150A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM150A CRISPRa kit - CRISPR gene activation of human transmembrane protein 150A
Symbol:	TMEM150A
Components:	GA115079G1, TMEM150A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115079G2, TMEM150A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115079G3, TMEM150A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q86TG1
Synonyms:	TM6PI; TMEM150
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

