


Product datasheet for **GA115080**

Human Nup53 (NUP35) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NUP35 CRISPRa kit - CRISPR gene activation of human nucleoporin 35
Symbol:	NUP35
Components:	GA115080G1, Nup53 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115080G2, Nup53 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115080G3, Nup53 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NFH5</u>
Synonyms:	MP44; NP44
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)