

Product datasheet for **GA115090**

Human TAF8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF8 CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 8
Symbol:	TAF8
Components:	GA115090G1, TAF8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115090G2, TAF8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115090G3, TAF8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z7C8</u>
Synonyms:	II; TAF; TAFII-43; TAFII43; TBN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)