

Product datasheet for **GA115091**

Human TMEM18 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM18 CRISPRa kit - CRISPR gene activation of human transmembrane protein 18
Symbol:	TMEM18
Components:	GA115091G1, TMEM18 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115091G2, TMEM18 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115091G3, TMEM18 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96B42
Synonyms:	DKFZp434C1714
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)