

Product datasheet for **GA115096**

Human TRIM43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRIM43 CRISPRa kit - CRISPR gene activation of human tripartite motif containing 43
Symbol:	TRIM43
Components:	GA115096G1, TRIM43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115096G2, TRIM43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115096G3, TRIM43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96BQ3
Synonyms:	TRIM43A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)