

## Product datasheet for **GA115104**

### Human REG3G activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	REG3G CRISPRa kit - CRISPR gene activation of human regenerating family member 3 gamma
Symbol:	REG3G
UniProt ID:	<a href="#">Q6UW15</a>
Synonyms:	LPPM429; PAP-1B; PAP1B; PAP IB; PAPIB; REG-III; REG III; UNQ429
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)