


Product datasheet for **GA115106**

Human CLHC1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CLHC1 CRISPRa kit - CRISPR gene activation of human clathrin heavy chain linker domain containing 1
Symbol:	CLHC1
Components:	GA115106G1, CLHC1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115106G2, CLHC1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115106G3, CLHC1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NHS4</u>
Synonyms:	C2orf63
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)