

Product datasheet for **GA115108**

Human AP1S3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AP1S3 CRISPRa kit - CRISPR gene activation of human adaptor related protein complex 1 subunit sigma 3
Symbol:	AP1S3
UniProt ID:	Q96PC3
Synonyms:	PSORS15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)