


Product datasheet for GA115115

Human KCTD18 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCTD18 CRISPRa kit - CRISPR gene activation of human potassium channel tetramerization domain containing 18
Symbol:	KCTD18
UniProt ID:	Q6PI47
Synonyms:	6530404F10Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)