

Product datasheet for **GA115124**

Human TMEM178 (TMEM178A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM178A CRISPRa kit - CRISPR gene activation of human transmembrane protein 178A
Symbol:	TMEM178A
UniProt ID:	<u>Q8NBL3</u>
Synonyms:	TMEM178
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)