

## Product datasheet for **GA115125**

### Human CPO activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                      |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                         |
| Description:  | CPO CRISPRa kit - CRISPR gene activation of human carboxypeptidase O                                                                                                                                                                                                                                 |
| Symbol:       | CPO                                                                                                                                                                                                                                                                                                  |
| Components:   | GA115125G1, CPO gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115125G2, CPO gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115125G3, CPO gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8IVL8</a></u>                                                                                                                                                                                                                                                                        |
| Synonyms:     | MGC138281; MGC138283                                                                                                                                                                                                                                                                                 |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                              |

#### Product images:



[View online »](#)