

Product datasheet for **GA115128**

Human **SLC66A3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC66A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 66 member 3
Symbol:	SLC66A3
Components:	GA115128G1, SLC66A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115128G2, SLC66A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115128G3, SLC66A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N755</u>
Synonyms:	C2orf22
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)