

## Product datasheet for **GA115136**

### Human **KCNH8** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNH8 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily H member 8
Symbol:	KCNH8
UniProt ID:	<a href="#">Q96L42</a>
Synonyms:	ELK; ELK1; elk3; Kv12.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)