

Product datasheet for **GA115151**

Human FAM43A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FAM43A CRISPRa kit - CRISPR gene activation of human family with sequence similarity 43 member A
Symbol:	FAM43A
Components:	GA115151G1, FAM43A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115151G2, FAM43A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115151G3, FAM43A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N2R8</u>
Synonyms:	FLJ90022
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)