

Product datasheet for **GA115153**

Human TMEM42 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM42 CRISPRa kit - CRISPR gene activation of human transmembrane protein 42
Symbol:	TMEM42
Components:	GA115153G1, TMEM42 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115153G2, TMEM42 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115153G3, TMEM42 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q69YG0</u>
Synonyms:	MGC29956
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)