

Product datasheet for **GA115159**

Human TMEM207 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM207 CRISPRa kit - CRISPR gene activation of human transmembrane protein 207
Symbol:	TMEM207
Components:	GA115159G1, TMEM207 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115159G2, TMEM207 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115159G3, TMEM207 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UWW9</u>
Synonyms:	UNQ846
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)