

Product datasheet for **GA115168**

Human SNTN activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SNTN CRISPRa kit - CRISPR gene activation of human sentan, cilia apical structure protein
Symbol:	SNTN
Components:	GA115168G1, SNTN gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115168G2, SNTN gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115168G3, SNTN gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A6NMZ2</u>
Synonyms:	S100A1L; S100AL; sentan
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)