

Product datasheet for **GA115174**

Human C4orf33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C4orf33 CRISPRa kit - CRISPR gene activation of human chromosome 4 open reading frame 33
Symbol:	C4orf33
UniProt ID:	Q8N1A6
Synonyms:	FLJ33703
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)