

Product datasheet for **GA115175**

Human SMIM43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM155 CRISPRa kit - CRISPR gene activation of human transmembrane protein 155
Symbol:	SMIM43
Components:	GA115175G1, SMIM43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115175G2, SMIM43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115175G3, SMIM43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4W5P6</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

