


Product datasheet for **GA115177**

Human TENR (ADAD1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ADAD1 CRISPRa kit - CRISPR gene activation of human adenosine deaminase domain containing 1
Symbol:	ADAD1
Components:	GA115177G1, TENR gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115177G2, TENR gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115177G3, TENR gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96M93</u>
Synonyms:	Tenr
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)