


Product datasheet for **GA115182**

Human **TMPRSS11B** activation kit by **CRISPRa**

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TMPRSS11B CRISPRa kit - CRISPR gene activation of human transmembrane serine protease 11B |
| Symbol: | TMPRSS11B |
| UniProt ID: | Q86T26 |
| Synonyms: | DKFZp686L1818 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)