

Product datasheet for **GA115222**

Human TMEM174 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TMEM174 CRISPRa kit - CRISPR gene activation of human transmembrane protein 174 |
| Symbol: | TMEM174 |
| UniProt ID: | Q8WUU8 |
| Synonyms: | FLJ31268; MGC13034 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)