

Product datasheet for **GA115251**

Human HUS1B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HUS1B CRISPRa kit – CRISPR gene activation of human HUS1 checkpoint clamp component B
Symbol:	HUS1B
Components:	GA115251G1, HUS1B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115251G2, HUS1B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115251G3, HUS1B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NHY5</u>
Synonyms:	RP11-532F6.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

