

Product datasheet for **GA115254**

Human WBS28 (WBSCR28) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM270 CRISPRa kit - CRISPR gene activation of human transmembrane protein 270
Symbol:	TMEM270
Components:	GA115254G1, WBS28 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115254G2, WBS28 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115254G3, WBS28 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6UE05
Synonyms:	MGC26719
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)