


Product datasheet for **GA115265**

Human LOC136242 (PRSS37) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRSS37 CRISPRa kit - CRISPR gene activation of human serine protease 37
Symbol:	PRSS37
Components:	GA115265G1, LOC136242 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115265G2, LOC136242 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115265G3, LOC136242 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A4D1T9</u>
Synonyms:	TRYX2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)