

Product datasheet for **GA115296**

Human DCAF4L2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DCAF4L2 CRISPRa kit - CRISPR gene activation of human DDB1 and CUL4 associated factor 4 like 2
Symbol:	DCAF4L2
Components:	GA115296G1, DCAF4L2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115296G2, DCAF4L2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115296G3, DCAF4L2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NA75</u>
Synonyms:	WDR21C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)