

## Product datasheet for **GA115310**

### Human TAF1L activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TAF1L CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 1 like
Symbol:	TAF1L
UniProt ID:	<a href="#">Q8IZX4</a>
Synonyms:	TAF(II)210; TAF2A2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

