


Product datasheet for **GA115330**

Human DGKK activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DGKK CRISPRa kit - CRISPR gene activation of human diacylglycerol kinase kappa
Symbol:	DGKK
Components:	GA115330G1, DGKK gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115330G2, DGKK gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115330G3, DGKK gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5KSL6</u>
Synonyms:	MRX68
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)