

Product datasheet for **GA115332**

Human PWWP3B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PWWP3B CRISPRa kit – CRISPR gene activation of human PWWP domain containing 3B
Symbol:	PWWP3B
Components:	GA115332G1, PWWP3B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115332G2, PWWP3B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115332G3, PWWP3B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5H9M0</u>
Synonyms:	FLJ33516; MGC129994; MGC129995
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

