


Product datasheet for **GA115348**

Human **GAB3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	GAB3 CRISPRa kit - CRISPR gene activation of human GRB2 associated binding protein 3
Symbol:	GAB3
Components:	GA115348G1, GAB3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115348G2, GAB3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115348G3, GAB3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WWW8</u>
Synonyms:	DOS/Gab family member 3; Gab3 protein; Gab3 scaffolding protein; GRB2-associated binding protein 3; OTTHUMP00000026055
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)