

Product datasheet for **GA115352**

Human GPR119 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | GPR119 CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor 119 |
| Symbol: | GPR119 |
| UniProt ID: | <u>Q8TDV5</u> |
| Synonyms: | GPCR2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)