

## Product datasheet for **GA115359**

### Human TCP10L activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TCP10L CRISPRa kit - CRISPR gene activation of human t-complex 10 like
Symbol:	TCP10L
UniProt ID:	<u><a href="#">Q8TDR4</a></u>
Synonyms:	C21orf77; PRED77; TCP10A-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)